



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Treatment table Wesseling BASIX 4 (BAS4EMW) (New)

New,
Dutch manufacture,
This two-piece treatment table is suitable for a variety of purposes.
BasiX 4 has height-adjustable armrests next to the head section,
Thanks to the armrests, the patient can put his hands next to his head and assume a very comfortable position,
As the head section is also negatively adjustable, the neck and shoulders can also be well positioned,
BasiX 4 is also available with a tapered leg section.
Thanks to the scissor frame, the table has a wide range of height adjustment and lifts straight up.
The headrest can be adjusted both positively and negatively.
As standard, the table is adjustable by means of an electric actuator and has wheels with braking system,
Technical parameters:
Length: 200 cm,
Width: 64 cm,
Upholstery thickness: 50 mm
Height adjustment: from 45 to 115 cm,
Headrest adjustment: 0° to 85°,
Adjustment of armrests: 0-17 cm,
Load capacity : 225 kg,
Safe working load: 250 kg,
Duration of height adjustment: 40 sec.
Adjustable sections: 3,
Brake system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Paper roll holder,
Nose hole,
Large neck cushion,
Foot height adjustment (control from the frame),
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
Terry cloth cover for tables 65-70 cm,
Terry cloth cover for tables 65-70cm with cut-out for face.
It has an up-to-date Technical Passport,
Warranty: 24 months

