



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Treatment table Wesseling BASIX 4T (BAS4TEMW) (New)**

New,  
Dutch manufacture,  
This two-piece treatment table is suitable for a variety of purposes.  
BasiX 4T has height-adjustable armrests next to the head section,  
Thanks to the armrests, the patient can put his hands next to his head and assume a very comfortable position,  
As the head section is also negatively adjustable, the neck and shoulders can also be well positioned,  
BasiX 4T is also available with a tapered leg section.  
Thanks to the scissor frame, the table has a wide range of height adjustment and lifts straight up.  
The headrest can be adjusted both positively and negatively.  
As standard, the table is adjustable by means of a hydraulic cylinder, has wheels with braking system and Trendelenburg position,  
Technical parameters:  
Length: 200 cm,  
Width: 64 cm,  
Upholstery thickness: 50 mm  
Height adjustment: 48 to 98 cm,  
Headrest adjustment: 0° to 85°,  
Sleep surface adjustment: 0° to -25°  
Adjustment of armrests: 0-17 cm,  
Load capacity : 225 kg,  
Safe working load: 250 kg,  
Duration of height adjustment: 40 sec.  
Adjustable sections: 4,  
Brake system: central,  
Wheels: Ø 75 mm,  
As additional options available (at extra charge):  
Paper roll holder,  
Nose hole,  
Large neck cushion,  
Foot height adjustment (control from the frame),  
Special table top width 70cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
Terry cloth cover for tables 65-70 cm,  
Terry cloth cover for tables 65-70cm with cut-out for face.  
It has an up-to-date Technical Passport,  
Warranty: 24 months

