



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling BASIX 5 (BAS5EMW) (New)

New,
Dutch production,
BASIX 5 is a 5-part table which combines all functions in one table.
The back and middle sections are equipped with two gas springs each, the headrest section with a single gas spring.
The table has height-adjustable armrests next to the head section
Thanks to the armrests, the patient can put his/her hands near the head and adopt a very comfortable position,
As standard, the table is adjustable with an electric actuator and has wheels with a braking system,
Technical parameters:
Length: 200 cm,
Width: 64 cm,
Upholstery thickness: 50 mm
Height adjustment: from 45 to 115 cm,
Headrest adjustment: -55° to 85°,
Backrest adjustment: -25° to 85°,
Adjustment of armrests: 0-17 cm,
Load capacity : 225 kg,
Safe working load: 250 kg,
Duration of height adjustment: 40 sec.
Adjustable sections: 4,
Brake system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Paper roll holder,
Nose hole,
Large neck cushion,
Foot height adjustment (control from the frame),
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
Terry cloth cover for tables 65-70 cm,
Terry cloth cover for tables 65-70cm with cut-out for face.
Additional switch for height adjustment,
It has a valid Technical Passport,
Warranty: 24 months

