



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Wesseling treatment table BOBATH 2 ELITE (BOB2BE110, BOB2BE120) (New)**

New,  
Dutch production,  
The Bobath 2 Elite is the largest two-section table in the Bobath range,  
The backrest section is fitted with two gas springs, which are operated by a single lever under the section itself.  
The tables are available in sizes 90 cm, 100 cm,  
Thanks to the new scissor mechanism, this model can reach the lowest position of 42 cm and the highest position of 112 cm,  
Thanks to its powerful 8000 N motor, the table can lift 200 kg as standard.  
Sturdy steel construction with epoxy paint coating,  
Height adjustment of the table is done by an electric actuator,  
The table is only available with castors and central locking,  
Technical parameters:  
Length: 200 cm,  
Width: 110 or 120 cm,  
Height adjustment: from 42 to 112 cm,  
Load capacity 200 kg,  
Safe working load 250 kg  
Duration of height adjustment (sec.): 40  
Central braking system,  
Wheels 75 Ø in mm,  
Adjustable sections 1,  
Headrest adjustment 0° to +85°,  
As optional extras available (at extra cost):  
Paper holder,  
Thicker upholstery 70 mm instead of 50 mm,  
Face hole in the headrest,  
It has a valid Technical Passport,  
24-month warranty

