



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling treatment table BOBATH 2 ELITE (BOB2BH110, BOB2BH120) (New)

New,
Dutch production,
The Bobath 2 Elite is the largest two-section table in the Bobath range,
The backrest section is fitted with two gas springs, which are operated by a single lever under the section itself.
The tables are available in sizes 110 cm, 120 cm,
Thanks to the new scissor mechanism, this model can reach the lowest position of 42 cm and the highest position of 112 cm,
Thanks to its powerful 8000 N motor, the table can lift 200 kg as standard.
Sturdy steel construction with epoxy paint coating,
Height adjustment of the table is done by a hydraulic cylinder,
The table is only available with castors and central locking,
Technical parameters:
Length: 200 cm,
Width: 90 or 100 cm,
Height adjustment: from 42 to 112 cm,
Load capacity 200 kg,
Safe working load 250 kg
Duration of height adjustment (sec.): 40
Central braking system,
Wheels 75 Ø in mm,
Adjustable sections 1,
Headrest adjustment 0° to +85°,
As optional extras available (at extra cost):
Paper holder,
Thicker upholstery 70 mm instead of 50 mm,
Face hole in the headrest,
It has a valid Technical Passport,
24-month warranty

