



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling ECO-2 ECG cardiology table (ECO3CGHMW) (New)

New,
Dutch production,
Thanks to the double hinges and the wide lower frame this table is very stable,
In the headrest there is a 30x30 cm cut-out with a plug which makes this table ideal for echocardiography,
Under the frame of the table there is a lot of space, which facilitates cooperation with the patient lift,
Sturdy steel structure coated with epoxy paint,
Height adjustment of the table is carried out by a hydraulic cylinder,
The table is only available with castors and a brake,
A two-part sleep surface is included as standard,
Technical parameters:
Length: 196 cm,
Width: 64 cm,
Height adjustment: from 46 to 96 cm,
Load capacity 175 kg,
Safe working load 225 kg
Duration of height adjustment: 26 sec.
Adjustable sections: 1,
Head restraint adjustment: 0° to 85°,
Central braking system,
Wheels 75 Ø in mm,
As additional options available (at extra charge):
Paper holder,
Thicker upholstery 70 mm instead of 50 mm,
Side rails,
Large neck cushion,
Foot height adjustment (control from the frame),
Special table top width 70cm,
Special top width 80cm,
It has an up-to-date Technical Passport,
Warranty: 24 months

