



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling folding 3-segment couch (KOF0003) (New)

New,
Dutch production,
Folding couch,
In standard there is a three-section couch bed made of plywood,
For an additional fee a cover for the couch,
Technical parameters:
Length: 197 cm (after folding 98 cm),
Width: 56 cm,
Fixed height: 78 cm,
Headrest adjustment: 0° to 90°,
Weight: 18 kg,
It has a valid Technical Passport,
Warranty: 24 months

