



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **Treatment table Wesseling MONDIAL 2 (MON2EMW) (New)**

---

New,  
Dutch production,  
The Mondial 2 is mounted on a special frame so the sections can be adjusted to almost any position, This makes this table ideal for physiotherapists.  
The headrest can be adjusted both positively and negatively.  
As standard, the table is adjusted using an electric actuator and has wheels with a braking system,  
The table has Trendelenburg, anti-Trendelenburg and drainage/flex positions.  
This model is equipped with a face cut-out in the headrest as standard,  
Technical parameters:  
Length: 200 cm,  
Width: 64 cm,  
Upholstery thickness: 50 mm  
Height adjustment: from 46 to 98 cm,  
Headrest adjustment: -20° to 75°,  
Load capacity : 175 kg,  
Safe working load: 225 kg,  
Duration of height adjustment: 26 sec.  
Adjustable sections: 2,  
Brake system: central,  
Wheels: Ø 75 mm,  
As additional options available (at extra charge):  
Paper roll holder,  
Large neck cushion,  
Foot height adjustment (controlled from the frame),  
Special table top width 70cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
It has a current Technical Passport,  
Warranty: 24 months

