



Treatment table Wesseling MONDIAL 4 (MON4HMW) (New)

New,
Dutch production,
The Mondial 4 is mounted on a special frame so the sections can be adjusted to almost any position, This makes this table ideal for physiotherapists.
The headrest can be adjusted both positively and negatively.
As standard, the table is adjustable via a hydraulic cylinder and has wheels with a braking system,
The table has Trendelenburg, anti-Trendelenburg and drainage/flex positions.
This model is equipped with a face cut-out in the headrest as standard,
Technical parameters:
Length: 200 cm,
Width: 64 cm,
Upholstery thickness: 50 mm
Height adjustment: from 46 to 98 cm,
Headrest adjustment: -90° to 15°,
Backrest section adjustment: -20° to 75°,
Back/body section adjustment: -25° to 75°,
Adjustment of leg section: -90° to 15°,
Load capacity : 175 kg,
Safe working load: 225 kg,
Duration of height adjustment: 27 sec.
Adjustable sections: 4,
Braking system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Paper roll holder,
Large neck cushion,
Foot height adjustment (controlled from the frame),
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
It has a current Technical Passport,
Warranty: 24 months

