



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling plastering treatment table PLASTER E (PLASH01) (New)

New,
Of Dutch manufacture,
Height adjustment of the table is done by a hydraulic cylinder,
Stable steel construction covered with epoxy paint,
In standard there is a four-section sleep surface made of plywood,
The foot section is bipartite,
With the long upright backrest a comfortable sitting position can be created,
All movable sections are adjustable with a gas spring,
The use of the table can be extended with various options such as monkey bars, knee/leg support and stainless steel rails on each side,
Technical specifications:
Length:
205 cm with patient holder on frame,
185 cm without patient handle,
Width: 60 cm,
Thickness of upholstery: 50 mm
Height adjustment: from 65 to 115 cm,
Backrest adjustment: 0° to 85°,
Adjustment of leg sections: -90° to 45°,
Load capacity : 175 kg,
Safe working load: 200 kg,
Duration of height adjustment: 26 sec.
Adjustable sections: 3,
Wheels: Ø 75 mm,
Brake system: central,
The following are available as optional extras (at extra cost):
MONKEY POLE,
Paper roll holder,
Large neck cushion,
Foot height adjustment (controlled from the frame),
It has a valid Technical Passport,
Warranty: 24 months

