



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Plastering treatment table Wesseling PLASTER COMPACT (PLACOM00) (New)

New,  
Of Dutch manufacture,  
Height adjustment of the table is done by an electric actuator,  
Stable steel construction covered with epoxy paint,  
In standard there is a four-section sleep surface made of plywood,  
The foot section is bipartite,  
With the long upright backrest a comfortable sitting position can be created,  
All movable sections are adjustable with a gas spring,  
The use of the table can be expanded with various options such as a monkey bar, knee/leg support and stainless steel rails on each side,  
Thanks to the round legs and thicker profiles, the table not only looks beautiful but is also very strong and stable,  
Technical parameters:  
Length: 196 cm,  
Width: 64 cm,  
Upholstery thickness: 50 mm  
Height adjustment: from 42 to 112 cm,  
Backrest adjustment: 0° to 85°,  
Adjustment of leg sections: -85° to 35°,  
Load capacity : 200 kg,  
Safe working load: 250 kg,  
Duration of height adjustment: 40 sec.  
Adjustable sections: 3,  
Wheels: Ø 75 mm,  
Brake system: central,  
As optional extras available (at extra cost):  
Paper roll holder,  
Large neck cushion,  
Foot height adjustment (controlled from the frame),  
Special table top width 70cm,  
It has a current Technical Passport,  
Warranty: 24 months

