



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Treatment table Wesseling PRACTICAL 3 (PRA3E01) (New)**

New,  
Of Dutch manufacture,  
In standard there is a three-section sleep surface made of plywood,  
Stable steel construction coated with epoxy paint,  
The adjustment of the table is carried out with an electric actuator and has wheels with a braking system,  
The table has Trendelenburg, drainage and flex-position.  
The table is designed so that the therapist can put his feet under the table to approach the patient as close as possible,  
Technical parameters:  
Length: 196 cm,  
Width: 70 cm,  
Upholstery thickness: 50 mm  
Height adjustment: from 46 to 96 cm,  
Adjustment of middle section: 0° to 20°,  
Headrest adjustment: -20° to 85°,  
Adjustment of leg/back section: -20° to 80°,  
Load capacity : 175 kg,  
Safe working load: 225 kg,  
Duration of height adjustment: 26 sec.  
Adjustable sections: 3,  
Brake system: central,  
Wheels: Ø 75 mm,  
As optional extras available (at extra cost):  
Face opening,  
Paper roll holder,  
Large neck cushion,  
Anti-slip pads,  
Special table top width 60cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
It has a current Technical Passport,  
Warranty: 24 months

