



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling PRACTICAL 3 (PRA3H01) (New)

New,
Of Dutch manufacture,
In standard there is a three-section sleep surface made of plywood,
Stable steel construction covered with epoxy paint,
The adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,
The table has Trendelenburg, drainage and flex-position.
The table is designed so that the therapist can put his feet under the table to approach the patient as close as possible,
Technical parameters:
Length: 196 cm,
Width: 70 cm,
Upholstery thickness: 50 mm
Height adjustment: from 46 to 96 cm,
Adjustment of middle section: 0° to 20°,
Headrest adjustment: -20° to 85°,
Adjustment of leg/back section: -20° to 80°,
Load capacity : 175 kg,
Safe working load: 225 kg,
Duration of height adjustment: 26 sec.
Adjustable sections: 3,
Brake system: central,
Wheels: Ø 75 mm,
As optional extras available (at extra cost):
Face opening,
Paper roll holder,
Large neck cushion,
Anti-slip pads,
Special table top width 60cm,
Special top width 80cm,
Thicker upholstery 70mm,
It has a current Technical Passport,
Warranty: 24 months

