



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling PRACTICAL 4 (PRA4HV1) (New)

New,
Of Dutch manufacture,
Standard features a four-section plywood sleep surface with a narrower top in the foot section,
Stable steel construction coated with epoxy paint,
The adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,
The table has Trendelenburg position,
The table is designed so that the therapist can put his legs under the table to approach the patient as close as possible,
Technical parameters:
Length: 196 cm,
Width: 70 cm,
Thickness of upholstery: 50 mm
Height adjustment: from 46 to 96 cm,
Adjustment of the main section: 0° to 35°,
Headrest adjustment: -20° to 85°,
Adjustment of armrests: 0° to 13 cm,
Load capacity : 175 kg,
Safe working load: 225 kg,
Duration of height adjustment: 26 sec.
Adjustable sections: 4,
Brake system: central,
Wheels: Ø 75 mm,
As optional extras available (at extra cost):
Face opening,
Paper roll holder,
Large neck cushion,
Anti-slip pads,
Special table top width 60cm,
Special top width 80cm,
Thicker upholstery 70mm,
It has a current Technical Passport,
Warranty: 24 months

