



## Treatment table Wesseling PRACTICAL 5 (PRA5H01) (New)

New,  
Of Dutch manufacture,  
Standard features five-section plywood sleep surface,  
Stable steel construction covered with epoxy paint,  
The adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,  
The table has Trendelenburg position, drainage/flex-position,  
The table is designed so that the therapist can put his legs under the table to approach the patient as close as possible,  
Technical parameters:  
Length: 196 cm,  
Width: 70 cm,  
Upholstery thickness: 50 mm  
Height adjustment: from 46 to 96 cm,  
Adjustment of the main section: 0° to 35°,  
Headrest adjustment: -20° to 80°,  
Adjustment of armrests: 0° to 13 cm,  
Load capacity : 175 kg,  
Safe working load: 225 kg,  
Duration of height adjustment: 26 sec.  
Adjustable sections: 5,  
Brake system: central,  
Wheels: Ø 75 mm,  
As optional extras available (at extra cost):  
Face opening,  
Paper roll holder,  
Large neck cushion,  
Anti-slip pads,  
Special table top width 60cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
It has a current Technical Passport,  
Warranty: 24 months

