



Treatment table Wesseling PRACTICAL 5 (PRA5H01) (New)

New,
Of Dutch manufacture,
Standard features five-section plywood sleep surface,
Stable steel construction covered with epoxy paint,
The adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,
The table has Trendelenburg position, drainage/flex-position,
The table is designed so that the therapist can put his legs under the table to approach the patient as close as possible,
Technical parameters:
Length: 196 cm,
Width: 70 cm,
Upholstery thickness: 50 mm
Height adjustment: from 46 to 96 cm,
Adjustment of the main section: 0° to 35°,
Headrest adjustment: -20° to 80°,
Adjustment of armrests: 0° to 13 cm,
Load capacity : 175 kg,
Safe working load: 225 kg,
Duration of height adjustment: 26 sec.
Adjustable sections: 5,
Brake system: central,
Wheels: Ø 75 mm,
As optional extras available (at extra cost):
Face opening,
Paper roll holder,
Large neck cushion,
Anti-slip pads,
Special table top width 60cm,
Special top width 80cm,
Thicker upholstery 70mm,
It has a current Technical Passport,
Warranty: 24 months

