



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Treatment table Wesseling PRACTICAL 5 (PRA5H01) (New)

New,  
Of Dutch manufacture,  
Standard features five-section plywood sleep surface,  
Stable steel construction covered with epoxy paint,  
The adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,  
The table has Trendelenburg position, drainage/flex-position,  
The table is designed so that the therapist can put his legs under the table to approach the patient as close as possible,  
Technical parameters:

Length: 196 cm,

Width: 70 cm,

Upholstery thickness: 50 mm

Height adjustment: from 46 to 96 cm,

Adjustment of the main section: 0° to 35°,

Headrest adjustment: -20° to 80°,

Adjustment of armrests: 0° to 13 cm,

Load capacity : 175 kg,

Safe working load: 225 kg,

Duration of height adjustment: 26 sec.

Adjustable sections: 5,

Brake system: central,

Wheels: Ø 75 mm,

As optional extras available (at extra cost):

Face opening,

Paper roll holder,

Large neck cushion,

Anti-slip pads,

Special table top width 60cm,

Special top width 80cm,

Thicker upholstery 70mm,

It has a current Technical Passport,

Warranty: 24 months

