



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling PRAXIS 3 ELITE treatment table (PRAX3EMWH) (New)

New,

Dutch production,

In standard there is a three-section sleep surface made of plywood,

Stable steel construction covered with epoxy paint,

Adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,

Technical parameters:

Length: 196 cm,

Width: 64 cm,

Thickness of upholstery: 50 mm

Height adjustment: from 48 to 98 cm

Head/leg section adjustment: 35 to -85°,

Backrest adjustment: 0° to 85°,

Adjustable sections: 2,

Load capacity : 225 kg,

Safe working load: 250 kg,

Duration of height adjustment: 40 sec.

Brake system: central,

Wheels: Ø 75 mm,

As additional options available (at extra charge):

Face opening,

Paper roll holder,

Large neck cushion,

Foot height adjustment (controlled from the frame),

Anti-slip pads,

Special table top width 70cm,

Special top width 80cm,

Thicker upholstery 70mm,

Adapters for arms / feet / footrests,

Armrest,

Set of footrests,

Stirrups (set),

Foot rests,

Comes with a current Technical Passport,

Warranty: 24 months

