



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling KOMBI ELITE (KOMEMWH) (New)

New,
Dutch production,
In standard there is a three-section sleep surface made of plywood,
Stable steel construction covered with epoxy paint,
The table is adjusted by means of a hydraulic cylinder and has wheels with a braking system,
With the use of a notch in the table top in the lower segment of the table and the installation of a container made of non-rusting steel, it can be used in gynaecological surgeries,
Technical parameters:
Length: 198 cm,
Width: 64 cm,
Thickness of upholstery: 50 mm
Height adjustment: from 48 to 98 cm
Head section adjustment: -35° to 85°,
Center/back section adjustment: 0° to 85°,
Adjustable sections: 2,
Lifting capacity : 225 kg,
Safe working load: 250 kg,
Duration of height adjustment: 40 sec.
Brake system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Face opening,
Paper roll holder,
Large neck cushion,
Foot height adjustment (controlled from the frame),
Anti-slip pads,
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
Adapters for arms / feet / footrests,
Armrest,
Set of footrests,
Stirrups (set),
Foot rests,
Comes with a current Technical Passport,
Warranty: 24 months

