



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LIPOSUNIC ULIPO device for non-surgical body shaping and slimming (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 200 W,

Liposonic device is a non-invasive fat tissue reduction, using highly concentrated ultrasound (HIFU technology),

With Ulipo, patients can experience true slimming and skin smoothing - by an average of 2.5 cm in the waistline - in just one hour of treatment, Ulipo causes temperatures to rise above 56°C, effectively destroying fat tissue without damaging the tissues outside the treatment area.

Macrophages (connective tissue cells) are attracted to the fat tissue being treated, attaching free lipids and cellular debris. The destroyed cells are processed naturally by the body and gradually released over a period of 8-12 weeks,

Ulipo is also an ideal device for people who are almost the perfect size and shape - but who need help when diet and exercise are no longer working,

Benefits of the device:

Targeted and effective body modelling tailored to the patient,

Results visible in every patient,

Non-surgical, non-invasive,

Safe and effective for everyone,

Included:

0.8 cm (4 Mhz) handpiece - for hand, leg and thigh treatments,

Handpiece 1.3 cm (4 Mhz) - for abdomen, thighs and buttocks treatments,

Instruction manual in Polish,

Dimensions: 46 x 35 x 100 cm,

Weight of set: 37 kg,

It has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): instalments, leasing, loan

