



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Wesseling BOBATH 2 ELITE rehabilitation table (BOB2BH110, BOB2BH120) (New)**

New,  
Dutch manufacture,  
The Bobath 2 Elite is the largest two-section table in the Bobath range,  
The backrest section is fitted with two gas springs, which are operated by a single lever under the section itself.  
The tables are available in sizes 110 cm, 120 cm,  
Thanks to the new scissor mechanism, this model can reach the lowest position of 42 cm and the highest position of 112 cm,  
Thanks to its powerful 8000 N motor, the table can lift 200 kg as standard.  
Sturdy steel construction with epoxy paint coating,  
Height adjustment of the table is done by a hydraulic cylinder,  
The table is only available with castors and central locking,  
Technical parameters:  
Length: 200 cm,  
Width: 90 or 100 cm,  
Height adjustment: from 42 to 112 cm,  
Load capacity 200 kg,  
Safe working load 250 kg  
Duration of height adjustment (sec.): 40  
Central braking system,  
Wheels 75 Ø in mm,  
Adjustable sections 1,  
Headrest adjustment 0° to +85°,  
As optional extras available (at extra cost):  
Paper holder,  
Thicker upholstery 70 mm instead of 50 mm,  
Face hole in the headrest,  
It has a valid Technical Passport,  
24-month warranty

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

