



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Rehabilitation table Wesseling ECO-2 ECG (ECO3CGEMW) (New)

New,  
Dutch production,  
Thanks to the double hinges and the wide lower frame this table is very stable,  
In the headrest there is a 30x30 cm cut-out with a plug which makes this table ideal for echocardiography,  
Under the frame of the table there is a lot of space, which facilitates cooperation with the patient lift,  
Sturdy steel structure coated with epoxy paint,  
Height adjustment of the table is carried out with an electric actuator,  
The table is only available with castors and a brake,  
A two-part sleep surface is included as standard,  
Technical parameters:  
Length: 196 cm,  
Width: 64 cm,  
Height adjustment: from 46 to 96 cm,  
Load capacity 175 kg,  
Safe working load 225 kg  
Duration of height adjustment: 26 sec.  
Adjustable sections: 1,  
Head restraint adjustment: 0° to 85°,  
Central braking system,  
Wheels 75 Ø in mm,  
As additional options available (at extra charge):  
Paper holder,  
Thicker upholstery 70 mm instead of 50 mm,  
Side rails,  
Large neck cushion,  
Foot height adjustment (control from the frame),  
Special table top width 70cm,  
Special top width 80cm,  
It has an up-to-date Technical Passport,  
Warranty: 24 months,  
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

