



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Wesseling PRAXIS 3 ELITE rehabilitation table (PRAX3EMWH) (New)**

New,  
Of Dutch manufacture,  
In standard there is a three-section sleep surface made of plywood,  
Stable steel construction covered with epoxy paint,  
Adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,  
Technical parameters:

Length: 196 cm,  
Width: 64 cm,  
Thickness of upholstery: 50 mm  
Height adjustment: from 48 to 98 cm  
Head/leg section adjustment: 35 to -85°,  
Backrest adjustment: 0° to 85°,  
Adjustable sections: 2,  
Load capacity : 225 kg,  
Safe working load: 250 kg,  
Duration of height adjustment: 40 sec.

Brake system: central,  
Wheels: Ø 75 mm,  
As additional options available (at extra charge):  
Face opening,  
Paper roll holder,  
Large neck cushion,  
Foot height adjustment (controlled from the frame),  
Anti-slip pads,  
Special table top width 70cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
Adapters for arms / feet / footrests,

Armrest,  
Set of footrests,  
Stirrups (set),  
Foot rests,  
Comes with a current Technical Passport,  
Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

