



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table Wesseling KOMBI ELITE (KOMEMWE) (New)

New,
Of Dutch manufacture,
In standard there is a three-section sleep surface made of plywood,
Stable steel construction coated with epoxy paint,
The adjustment of the table is carried out with an electric actuator and has wheels with a braking system,
With the use of a notch in the table top in the lower segment of the table and the installation of a container made of non-rusting steel, it can be used in gynaecological surgeries,
Technical parameters:
Length: 198 cm,
Width: 64 cm,
Thickness of upholstery: 50 mm
Height adjustment: from 42 to 112 cm
Head section adjustment: -35° to 85°,
Center/back section adjustment: 0° to 85°,
Adjustable sections: 2,
Lifting capacity : 225 kg,
Safe working load: 250 kg,
Duration of height adjustment: 40 sec.
Brake system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Face opening,
Paper roll holder,
Large neck cushion,
Foot height adjustment (controlled from the frame),
Anti-slip pads,
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
Adapters for arms / feet / footrests,
Armrest,
Set of footrests,
Stirrups (set),
Foot rests,
Comes with a current Technical Passport,
Warranty: 24 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

