



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling rehabilitation table MONDIAL 3 (MON3EMW) (New)

New,
Dutch manufacture,
The Mondial 3 is mounted on a special frame so the sections can be adjusted to almost any position, This makes this table ideal for physiotherapists.
The headrest can be adjusted both positively and negatively.
As standard, the table is adjusted using an electric actuator and has wheels with a braking system,
The table has Trendelenburg, anti-Trendelenburg and drainage/flex positions.
This model is equipped with a face cut-out in the headrest as standard,
Technical parameters:
Length: 200 cm,
Width: 64 cm,
Upholstery thickness: 50 mm
Height adjustment: from 46 to 98 cm,
Headrest adjustment: -90° to 15°,
Backrest section adjustment: -20° to 75°,
Back/body section adjustment: -25° to 75°,
Load capacity : 175 kg,
Safe working load: 225 kg,
Duration of height adjustment: 26 sec.
Adjustable sections: 3,
Brake system: central,
Wheels: Ø 75 mm,
As additional options available (at extra charge):
Paper roll holder,
Large neck cushion,
Foot height adjustment (controlled from the frame),
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
It has a current Technical Passport,
Warranty: 24 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

