



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Wesseling rehabilitation table MONDIAL 4 (MON4HZW) (New)**

New,  
Dutch manufacture,  
The Mondial 4 is mounted on a special frame so the sections can be adjusted to almost any position, This makes this table ideal for physiotherapists.  
The headrest can be adjusted both positively and negatively.  
As standard, the table is adjusted using a hydraulic cylinder without wheels,  
The table has Trendelenburg, anti-Trendelenburg and drainage/flex positions.  
This model is equipped with a face cut-out in the headrest as standard,  
Technical parameters:  
Length: 200 cm,  
Width: 64 cm,  
Upholstery thickness: 50 mm  
Height adjustment: from 46 to 98 cm,  
Headrest adjustment: -90° to 15°,  
Backrest section adjustment: -20° to 75°,  
Back/body section adjustment: -25° to 75°,  
Adjustment of leg section: -90° to 15°,  
Load capacity : 175 kg,  
Safe working load: 225 kg,  
Duration of height adjustment: 27 sec.  
Adjustable sections: 4,  
Braking system: central,  
As additional options available (at extra charge):  
Paper roll holder,  
Large neck cushion,  
Foot height adjustment (control from the frame),  
Special table top width 70cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
It has a current Technical Passport,  
Warranty: 24 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

