



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling rehabilitation table MONDIAL 6 (MON6HMW) (New)

New,
Dutch manufacture,
The Mondial 6 is mounted on a special frame so the sections can be adjusted to almost any position, This makes this table ideal for physiotherapists.

The headrest can be adjusted both positively and negatively.

As standard, the table is adjustable via a hydraulic cylinder and has wheels with a braking system,

The table has Trendelenburg, anti-Trendelenburg and drainage/flex positions.

This model comes standard with adjustable armrests and a face cut-out in the headrest,

Technical parameters:

Length: 200 cm,

Width: 64 cm,

Upholstery thickness: 50 mm

Height adjustment: from 46 to 96 cm,

Headrest adjustment: -90° to 15°,

Back section adjustment: -20° to 75°,

Seat section adjustment: -25° to 75°,

Adjustment of armrests: 0° to -16°,

Load capacity : 175 kg,

Safe working load: 225 kg,

Duration of height adjustment: 27 sec.

Adjustable sections: 6,

Braking system: central,

Wheels: Ø 75 mm,

As additional options available (at extra charge):

Paper roll holder,

Large neck cushion,

Foot height adjustment (controlled from the frame),

Special table top width 70cm,

Special top width 80cm,

Thicker upholstery 70mm,

It has a current Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

