



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling table PLASTER MASTER (PLASTMAST) (New)

New,
Dutch manufacture,
The table features a sophisticated design with thoughtful details and modular accessories,
Height adjustment of the table is done by an electric actuator,
Sturdy steel construction with epoxy paint coating,
The table is characterised by a sophisticated design with well thought-out details and modular accessories,
Table advantages:

Easy access to the controls,
Very sturdy frame made of heavy metal components,
High stability - The table rests on four legs during treatment,
The table becomes mobile by means of a foot control,
Choice of wheel sizes 75 mm or 125 mm,
Table legs with anti-slip layer,

Technical parameters:

Length: 212cm with horizontal legs / 165cm with vertical legs,

Width: 64cm,

Upholstery thickness: 50 mm

Height adjustment: from 61 to 111cm,

Backrest adjustment: 0° to 90°,

Adjustment of thigh section: 0° to 65°,

Adjustment of shank section: 45° to -90°,

Load capacity : 175 kg,

Safe working load: 200 kg,

Adjustable sections: 3,

Wheels: Ø 75 mm,

Braking system: central,

Weight of equipment: 150 kg,

As additional options are available (for an extra charge):

MONKEY POLE,

Paper roll holder,

Large neck cushion,

Foot height adjustment (controlled from the frame),

It has a valid Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

