



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Wesseling table PLASTER E (PLASH01) (New)**

New,  
Of Dutch manufacture,  
Height adjustment of the table is done by a hydraulic cylinder,  
Stable steel construction covered with epoxy paint,  
In standard there is a four-section sleep surface made of plywood,  
The foot section is bipartite,  
With the long upright backrest a comfortable sitting position can be created,  
All movable sections are adjustable with a gas spring,  
The use of the table can be extended with various options such as monkey bars, knee/leg support and stainless steel rails on each side,  
Technical specifications:

Length:

205 cm with patient holder on frame,

185 cm without patient handle,

Width: 60 cm,

Thickness of upholstery: 50 mm

Height adjustment: from 65 to 115 cm,

Backrest adjustment: 0° to 85°,

Adjustment of leg sections: -90° to 45°,

Load capacity : 175 kg,

Safe working load: 200 kg,

Duration of height adjustment: 26 sec.

Adjustable sections: 3,

Wheels: Ø 75 mm,

Brake system: central,

The following are available as optional extras (at extra cost):

MONKEY POLE,

Paper roll holder,

Large neck cushion,

Foot height adjustment (controlled from the frame),

It has a valid Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

