



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Rehabilitation table Wesseling PRACTICAL 3 (PRA3E01) (New)

New,
 Dutch production,
 In standard there is a three-section sleep surface made of plywood,
 Stable steel construction coated with epoxy paint,
 The adjustment of the table is carried out with an electric actuator and has wheels with a braking system,
 The table has Trendelenburg, drainage and flex-position.
 The table is designed so that the therapist can put his feet under the table to approach the patient as close as possible,
 Technical parameters:

- Length: 196 cm,
- Width: 70 cm,
- Upholstery thickness: 50 mm
- Height adjustment: from 46 to 96 cm,
- Adjustment of middle section: 0° to 20°,
- Headrest adjustment: -20° to 85°,
- Adjustment of leg/back section: -20° to 80°,
- Load capacity : 175 kg,
- Safe working load: 225 kg,
- Duration of height adjustment: 26 sec.
- Adjustable sections: 3,
- Brake system: central,
- Wheels: Ø 75 mm,
- As optional extras available (at extra cost):
- Face opening,
- Paper roll holder,
- Large neck cushion,
- Foot height adjustment (controlled from the frame),
- Anti-slip pads,
- Special table top width 60cm,
- Special top width 80cm,
- Thicker upholstery 70mm,
- It has a current Technical Passport,
- Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE



WYBÓR BARW

