



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Rehabilitation table Wesseling PRACTICAL 4 (PRA4H01) (New)**

New,  
Dutch production,  
Standard features four-part plywood sleep surface,  
Stable steel construction covered with epoxy paint,  
The adjustment of the table is carried out with a hydraulic cylinder and has wheels with a braking system,  
The table has Trendelenburg position,  
The table is designed so that the therapist can put his legs under the table to approach the patient as close as possible,  
Technical parameters:  
Length: 196 cm,  
Width: 70 cm,  
Thickness of upholstery: 50 mm  
Height adjustment: from 46 to 96 cm,  
Adjustment of the main section: 0° to 35°,  
Headrest adjustment: -20° to 85°,  
Adjustment of armrests: 0° to 13 cm,  
Load capacity : 175 kg,  
Safe working load: 225 kg,  
Duration of height adjustment: 26 sec.  
Adjustable sections: 4,  
Brake system: central,  
Wheels: Ø 75 mm,  
As optional extras available (at extra cost):  
Face opening,  
Paper roll holder,  
Large neck cushion,  
Foot height adjustment (controlled from the frame),  
Anti-slip pads,  
Special table top width 60cm,  
Special top width 80cm,  
Thicker upholstery 70mm,  
It has a current Technical Passport,  
Warranty: 24 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE



PRAISTON



PRAISTON




PRAISTON



PRAISTON



PRAISTON



PRAISTON



PRAISTON