



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wesseling upright table ARISTO (ARIST00) (New)

New,
Dutch production,
Aristo is a very strong and rigid tilt table,
The table has two electric motors that enable adjustment of the couch to 90°,
Standard features include a two-part plywood sleep surface and wheels with central braking system,
Adjustment by means of electric actuators powered from a battery, its high capacity battery can be recharged by means of an adapter included in the set,
The patient can be secured to the table with optional straps, so the patient is safe and the table can be tilted to any angle from 0 ° to 90 °
What angle the table is set at can be seen on the protractor mounted on the side of the table,
The headrest is adjustable in any position by means of a gas spring. The footrests are individually adjustable in 6 positions. In position 6 the footrests are completely folded down,
Technical parameters:
Length: 220 cm,
Width: 85 cm,
Upholstery thickness: 50 mm
Height adjustment: from 53 to 103 cm,
Adjustment of the headrest: -35° to 90°,
Footrest adjustment: -20° to 20°,
Load capacity : 175 kg,
Safe working load : 175 kg,
Tilt adjustment duration: 60 sec.
Duration of height adjustment: 40 sec.
Adjustable sections: 1,
Brake system: central,
Wheels: Ø 100 mm,
As additional options are available (at extra cost):
Working platform for Aristo with height and depth adjustment,
Working platform for Aristo with height, depth and angle adjustment,
Handle set for Aristo adjustable in height,
Armrest for Aristo adjustable in height and depth, with adjustable handles,
Footrest for Aristo - one section,
Plastic leg bases for mounting on the Aristo footrest,
Fastening belt with Velcro for fixing the patient to the treatment table (180x8 cm),
Fastening belt with Velcro for fixing the patient to the treatment table (240x8 cm),
Cushion for the fastening belt provides additional comfort,
Linak battery for Jumbo - pack,
Paper roll holder,
Nose hole,
Large neck cushion,
Special top width 70cm,
Thicker upholstery 70mm,
Comes with a current Technical Passport issued,
Warranty: 24 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

