



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Endoscopic table Wesseling PRACTICAL 2 (PRA2ET0) (New)

New,
Dutch production,
In standard there is a two-part sleep surface made of plywood,
Both parts can be adjusted separately, in any position,
Sturdy steel structure coated with epoxy paint,
The adjustment of the table is carried out with an electric actuator and has wheels with a braking system,
The table has a Trendelenburg position,
Technical parameters:
Length: 196 cm,
Width: 70 cm,
Thickness of upholstery: 50 mm
Height adjustment: from 46 to 98 cm,
Headrest adjustment: -20° to 85°,
Adjustment of the main section: 0° to 35°,
Load capacity : 175 kg,
Safe working load: 225 kg,
Duration of height adjustment: 26 sec.
Adjustable sections: 2,
Brake system: central,
Wheels: Ø 75 mm,
As optional extras available (at extra cost):
Face opening,
Paper roll holder,
Large neck cushion,
Foot height adjustment (controlled from the frame),
Side rails,
Anti-slip pads,
Special table top width 70cm,
Special top width 80cm,
Thicker upholstery 70mm,
It has a current Technical Passport,
Warranty: 24 months,

