



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Endoscopic table Wesseling ECO-LINE 3 (ECO3EMW) (New)**

New,  
Dutch production,  
Thanks to double hinges and wide lower frame this table is very stable,  
There is a lot of space under the table frame which makes it easy to work with a patient lift,  
The table is equipped with an adjustable backrest section, thanks to which the Eco 3 can be placed in a sitting position,  
Sturdy steel structure coated with epoxy paint,  
Height adjustment of the table is carried out by means of an electric actuator,  
The table is only available with castors and a brake,  
In standard there is a three-part plywood recliner with a face hole,  
Technical parameters:  
Length: 196 cm,  
Width: 64 cm,  
Height adjustment: from 46 to 96 cm,  
Load capacity 175 kg,  
Safe working load 225 kg  
Duration of height adjustment: 26 sec.  
Adjustable sections: 2,  
Headrest adjustment: -30° to 85°,  
Backrest adjustment: 0° to 85°,  
Central braking system,  
Wheels 75 Ø in mm,  
As additional options available (at extra charge):  
Paper holder,  
Thicker upholstery 70 mm instead of 50 mm,  
Side rails,  
Large neck cushion,  
Foot height adjustment (control from the frame),  
Special table top width 70cm,  
Special top width 80cm,  
It has an up-to-date Technical Passport,  
Warranty: 24 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

