



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Endoscopic table Wesseling MONDIAL 5 (MON5EMW) (New)

New,  
Dutch production,  
Mondial 5 is mounted on a special frame so that the sections can be adjusted to almost any position, This makes this table ideal for physiotherapists.

The headrest can be adjusted both positively and negatively.

As standard, the table is adjusted using an electric actuator and has wheels with a braking system,

The table has Trendelenburg, anti-Trendelenburg and drainage/flex positions.

This model comes standard with adjustable armrests and a face cut-out in the headrest,

Technical parameters:

Length: 200 cm,

Width: 64 cm,

Upholstery thickness: 50 mm

Height adjustment: from 46 to 98 cm,

Headrest adjustment: -90° to 15°,

Adjustment of the middle section: -20° to 75°,

Adjustment of armrests: 0° to -16°,

Load capacity : 175 kg,

Safe working load: 225 kg,

Duration of height adjustment: 27 sec.

Adjustable sections: 5,

Brake system: central,

Wheels: Ø 75 mm,

As additional options available (at extra charge):

Paper roll holder,

Large neck cushion,

Foot height adjustment (controlled from the frame),

Special table top width 70cm,

Special top width 80cm,

Thicker upholstery 70mm,

It has a current Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

