



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**Autoclave BRAVO G4 28l (New)**

New  
Canadian production,  
Class B autoclave,  
Autoclave BRAVO G4, is the ultra-modern way to efficiently, quickly and reliably sterilise medical instruments in the office, clinic and hospital,  
The autoclave is managed entirely by an advanced on-board computer, with a spacious sterilization chamber made of pressed stainless steel,  
computer network support, cycle recording and other modern solutions placing it among the leading constructions of this type available on the market,  
It is characterised by a very advanced fractional vacuum system, which allows for complete (100%) air removal from hollow and porous materials, as well as an effective vacuum final drying phase. Removes all traces of moisture from any load,  
A powerful steam generation system, an extremely efficient hydraulic system and microprocessor control (combined with high-precision sensors), guarantee high process speed and excellent stability of thermodynamic parameters,  
Standard equipment BRAVO G4 28L:  
built-in cycle recorder,  
USB flash drive for the recorder,  
container for filling the clean water tank,  
rotating rack for tool trays,  
tool trays 5 pcs,  
funnel  
drain hose for drainage system,  
set of connectors,  
HEPA filter (installed in the device,  
operating manual in Polish,  
Technical data:  
Manufacturer: SCICAN,  
Standards: PN EN 13060,  
Application of the autoclave: medical autoclave,  
CE certificate: yes,  
Notification to UDT: full notification,  
Medical device certificate: yes,  
Warranty: 24 months or 2500 cycles,  
Language of use: Polish,  
Display: colour touchscreen,  
Registration of cycles: built-in recorder on pendrive,  
Sterilization of batches of type: B, S, N,  
Label printing: external printer (option),  
LAN operation: WiFi or RJ45 network cable,  
Network software: built-in, web browser,  
Chamber volume: 28l,  
Door locking: automatic,  
Chamber dimension: fi 28 cm length 45 cm,  
Number of trays: 6,  
Tray dimension: 18 x 38 cm,  
Maximum loading capacity:  
9 kg utensils,  
1,5 kg textiles  
Steam generator: outside the sterilisation chamber,  
Distilled water container: yes,  
Waste water tank: yes,  
External demineralizer: option,  
Connection to mains water and sewage system: option,  
Volume of the device: < 67dB,  
Power consumption: 2300W,  
External dimensions of the device: 48.0 x 60.0 x 45.0 cm (W x D x H),  
Weight of appliance: 53 kg,  
Emergency unlocking and opening of the appliance: yes,  
Additional technologies:  
LAN monitoring and management,  
G4 CLOUD - cloud operation,  
pre-programmed start-up,  
user profiles,  
test reminder function (Vacuum - Helix / B&D - Vacuum + Helix / B&D),  
It has an up-to-date technical passport,  
Warranty: 24 months or 2500 cycles (Warranty expires 24 months after purchase or after 2500 sterilization cycles)  
Possibility of financing (Poland only): instalments, leasing, loan.  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

