



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medical laser QUANTA SYSTEM Q1 LIGHT B (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Italian production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 3500 W,
- **The LIGHT B laser** is designed not only for a full range of vascular treatments, but also for hair removal and skin rejuvenation,
- **The device is equipped with:**
  - Built-in 5.5" touch screen,
  - Simple, intuitive menu, language selection: EN, IT,
  - Knob for changing parameter values,
  - **Laser source emitting a wavelength of 1064 nm:**
    - The 1064 nm Nd:Yag laser is a high-power laser with short and long pulse modes for effective vascular treatment, hair removal for all skin types and skin rejuvenation.
    - Due to its high absorption in hemoglobin, the 1064 nm wavelength is ideal for vascular treatment and discoloration removal,
    - The 1064 nm wavelength also has good absorption in melanin, which guarantees effectiveness in hair removal,
    - The short pulse mode is responsible for skin rejuvenation,
  - Safety button "STOP",
  - Fast TWAIN IPL calibration port,
  - Fiber optic,
  - Skin cooling tip (OFF, + 4° C, + 8° C, + 12° C),
- **Application:**
  - Permanent hair reduction for all skin phototypes,
  - Dermatological vascular treatments,
  - Skin rejuvenation treatments,
- **Laser specifications:**
  - Laser: ND:YAG,
  - Wavelength (nm):
    - 1064 long pulse,
    - 1064 short pulse,
  - Pulse length:
    - 3 - 300 ms (long pulse),
    - 0.3 - 0.8 ms (short pulse),
  - Repetition rate:
    - Up to 2 Hz (long pulse),
    - Up to 7 Hz (short pulse),
  - Spot size:
    - Up to 14 mm (long pulse),
    - Up to 6 mm (short pulse),
  - Maximum fluency [J/cm<sup>2</sup>]:
    - >800 @ 2 mm (long pulse),
    - Up to 30 @ 5 mm (short pulse),
- The mobile camera has two double, swivel wheels at the front,
- Convenient handles on the front and back for moving the device,
- **Included in the set:**
  - Laser head with a skin cooling tip,
  - New optical fiber,
  - Foot attachment,
  - Device keys - 2 pcs.,
  - Head holder,
  - User manual in Polish and English,
  - Power cord,
- Overall dimensions: 48 x 156 x 109 cm,
- Set weight: 130 kg,
- Has a valid Technical Passport issued,
- Guarantee:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Financing options (Poland only): Installments, Leasing, Loan,

