



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**Medical laser QUANTA SYSTEM Q1 LIGHT B (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Made in Italy,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 3500 W,  
Laser LIGHT B is designed not only for a full range of vascular treatments, but also for depilation and skin rejuvenation,  
The device is equipped with:  
Built-in 5.5 " touch screen,  
Simple intuitive menu, choice of language: EN, IT,  
Knob for changing parameter values,  
Laser source, emitting 1064 nm wavelength:  
1064 nm Nd:Yag laser, is a high power laser with short and long pulse modes for effective vascular treatment, hair removal from all skin types and skin rejuvenation.  
Due to its high absorption in haemoglobin, the 1064 nm wavelength is ideal for vascular treatment and pigmentation removal,  
The 1064 nm wavelength also has good absorption in melanin, which guarantees effectiveness in hair removal,  
Short pulse mode is responsible for skin rejuvenation,  
Safety "STOP" button,  
High speed TWAIN IPL calibration port,  
Light guide,  
Skin cooling tip (OFF, + 4° C, + 8° C, + 12° C),  
Application:  
Permanent hair reduction of any skin phototype,  
Dermatological vascular treatments,  
Skin rejuvenation treatments,  
Laser specifications:  
Laser: ND:YAG,  
Wavelength (nm):  
1064 long pulse,  
1064 short pulse,  
Pulse length:  
3 - 300 ms (long pulse),  
0.3 - 0.8 ms (short pulse),  
Repetition rate:  
Up to 2 Hz (long pulse),  
Up to 7 Hz (short pulse),  
Spot size:  
Up to 14 mm (long pulse),  
Up to 6 mm (short pulse),  
Maximum fluidity [J/cm2]:  
>800 @ 2 mm (long pulse),  
Up to 30 @ 5 mm (short pulse),  
Mobile apparatus, has two double swivel castors at the front,  
Convenient handles at the front and back to move the unit,  
Included:  
Laser head with skin cooling tip,  
Light guide,  
Foot attachment,  
Keys for the device - 2 pieces,  
Handle for the head,  
User manual in Polish and English,  
Power cord,  
Overall dimensions: 48 x 156 x 109 cm,  
Weight of set: 130 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): instalments, leasing, loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

