



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HP charger for CodeMaster 100 M2480B defibrillator batteries (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Manufacture U.S.A,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 200 VA,

The HP2480B charger is used to charge M2476 and M2477B batteries for the HP CodeMaster 100 defibrillator,

Allows two batteries (not supplied) to be charged simultaneously,

Equipped with a male RS-232 connector,

Features lights for charging and recharging batteries, performance test in progress, successful test, and error,

Keypad on control panel in French,

A visual alarm indicates when a battery test is due,

A legend that slides out from under the unit for easy operation,

Specifications:

Battery charging time: 2 - 3 h,

Performance test time: approx. 24 h,

Recommended charging temperature: 10 - 30° C,

Included:

English user manual (PDF),

Power cord,

Dimensions: 36.5 x 13.5 x 34 cm,

Weight: 7 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

