



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Shower stool ETAC STAPEL (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Swedish production,

ETAC Stapel shower stool is an ideal product for elderly people, making everyday hygiene easier,

The stool ensures safety thanks to rubber endings on the legs,

Anti-corrosive frame made of powder-coated steel,

Light, stable and durable construction,

The seat surface is equipped with a special structured layer which has an anti-slip effect,

Removable seat,

Mechanically washable at maximum 85° C for 3 minutes,

Dimensions:

Seat: 43 x 38 cm,

Seat height from the floor: 61 cm,

Overall dimensions: 49 x 61 x 49 cm,

Maximum load capacity: 130 kg,

Weight: 4.2 kg,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE



