



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Light source OLYMPUS CLV- S20 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 300 W,
OLYMPUS CLV - S20 is a xenon light source used in diagnostic and surgical procedures,
Manual and automatic mode,
Xenon light,
Bulbs:
Xenon 300 W Model MD-631 with short arc ozone-free - main lighting,
Halogen 150 W - emergency lighting,
Lamp life:
Xenon: 300 h of continuous use,
Halogen: 500 h,
Bulb life indicator,
18-stage regulation of light intensity,
Foiled control panel,
Easy access to the bulbs on the side of the unit,
Colour temperature: 6000 K,
Has inputs and outputs:
VIDEO IN/OUT,
25-pin female connector SYSTEM,
LIGHT CONTROL (cable not included, can be retrofitted at additional cost),
Includes:
Power cord,
Dimensions: 37 x 14 x 42 cm,
Weight: 12.1 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

