



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDI CAPTURE MVR Pro Medical Recorder (New)

New,
U.S.A. production,
Features and benefits of the Medi Capture MVR Pro recorder:
The MVR Pro HD is a highly innovative, yet easy-to-use, full HD medical video recorder,
It includes a spacious built-in touchscreen monitor and advanced networking features that are essential for file sharing in doctors' offices, health systems and their facilities,
Records Full HD video at 60 frames per second with smooth, crystal-clear resolution,
Captures still frames during recording,
Pre-recording function ensures that important recorded procedures are not missed,
Plays back in real time to check video (but not for diagnostic purposes),
Multi-input recording provides procedures for two cameras with easy switching function,
Pre-recording functionality ensures users don't miss important events, and the review function allows users to quickly check the recording for an extra level of security,
This highly advanced medical recorder is intuitive and requires no formal training,
Its lightweight, compact size easily connects to a camera and HD video monitor and can be operated by:
AMOLED 7"/18 cm integrated touchscreen monitor,
MediCapture's user-friendly icons,
Front panel buttons,
Remote control (included),
Pedal,
The MVR Pro HD is compact, easy to install and rugged enough for any operating room situation with its stainless steel housing and fluid resistant front panel,
It automatically prompts you to enter patient information via the keyboard into the monitor before recording. Patient information is saved along with the recording to a USB or external hard drive
for efficient workflow and data security,
Videos and images are saved along with patient data on a USB flash drive (or external USB hard drive), internal and network memory for archiving or viewing after surgery,
Kit contents:
MVR Pro recorder,
USB flash drive,
HDMI cable,
Wi-Fi / BT antenna,
HDMI to DVI adapter,
Adapter for split remote control cable,
Power cable,
Rubber feet (4),
User manual (English),
Additional accessories (extra charge):
DICOM activation key for DICOM Worklist/DICOM Store,
4K activation key for 2160p30 recording,
White slimcool+ keyboard with touchpad, English layout,
64GB USB flash drive,
512GB internal storage capacity instead of 128GB,
Internal storage capacity of 1TB instead of 128GB,
Technical specifications:
Device type: HD medical video recorder. Records high definition video and still images from medical video sources to external USB drive, network storage, Exports video and images to PACS in DICOM format,
Housing / Buttons: Stainless steel, Fluid resistant,
User interface: Internal touch screen overlay or external monitor
Video inputs:
2 x HDMI, DVI (with adapter),
4K-SDI with MVC Pro SDI to HDMI (optional),
Video outputs: 1x HDMI, DVI (with adapter),
Video input resolutions:
4K input only: 3840x2160p30/25,
Both inputs (1+2): 1920x1200p60
1080p60/50/30/25,
1080i60/50,
1280x1024p60/50,
720p60/50,
Video recording resolutions:
3840x2160,
1920x1200,
1920x1080,
1280x1024,
1280x720,
*For 3840x2160 recording, activate "4K recording" (option),
Video recording speed: MP4/coding H.264 up to 32 Mbit/s (1920x1080 p60, lower resolutions have lower bitrate), MP4/coding HEVC (H.265) reduces bitrate by 50%,
Photo format: JPG,PNG,
Audio Recording: Bluetooth microphone, stereo microphone (4 pole), USB microphone / sound embedded in video,
Internal memory: 128 GB. Expansion on request (at extra charge)
External drive: USB flash drive, USB hard drive (FAT32, NTFS, exFAT),
USB support: USB 3.0, super speed:
1 front connector,
2 rear connectors,
Network Storage: External server, 1000/100/10Mbit/s,
Footswitch/Trigger: 3.5mm mini jack socket. Compatible with standard foot switches (not included) and closed contact camera triggers, Bluetooth foot switch.
Remote control: infrared remote control (included),
Size (closed touch screen):
Size: 173x153x50mm/ 6.7x6.0x2.0",
Weight: 1.2kg,
Work Environment:
Storage: -40 ° to + 85 ° C,
Operation: -20 ° to +40 ° C,
Power supply:
AC 100-240V 50/60Hz 0.4-0.2A,
3-pin AC socket,
Certifications:
Class 1 equipment according to EU MDR 2017/745 and US-FDA,
EN 60601-1:2006 + A1:2013 + A12:2014 (Edition 3.1),
EMC: IEC 60601-1-2:2014 (4th edition), EN 60601-1-2:2015,
Utility: EN 60601-1-6:2010+A1:2015,
Radio equipment: RED 2014/53/EU,
FCC part 15,
WEEE / RoHS3 / ISO 13485: 2016
It has a current Technical Passport issued,
Warranty 24 months,
Financing option (Poland only): instalments, leasing, loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

