



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient monitor CRITICARE COMFORTCUFF 506N3 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
U.S.A. production,
Power supply: 230 V,
Emergency power supply: 1 built-in 6 V 7.2 Ah battery (not supplied) - can be retrofitted at additional cost
Frequency: 50/60 Hz,
Power: 60 VA,
The Criticare ComfortCuff is a lightweight, compact and mobile patient monitor ideal for spot measurements or continuous observation of vital functions,
Product features:
ComfortCuff technology of measuring pressure as the cuff is inflated allows for a more accurate result in less time,
Displays:
3.5" text LCD,
Digital LEDs: for systolic and diastolic blood pressure, mean blood pressure, saturation and pulse,
10 LED heart rate bar provides a graphic expression of the heart rate,
Clear and easy to use foiled control panel,
English keyboard on the control panel,
Menu selection in EN, DE, FR, IT, PT, ES, DA, NL,
Customisable audible and visual alarms,
10 levels of intensity for the audible alarm,
Option to silence the alarm for 2 minutes or permanently,
Indicator lights for low battery and wired power supply,
Possibility of storing 24 hours of measurements in the device memory,
Comfortable handle for carrying the monitor,
Specifications:
NIBP:
Measurement range:
Adults: 30 - 300 mmHg,
Children: 30 - 150 mmHg,
Newborns: 20 - 150 mmHg,
Measurement cycle: < 40 s,
Measurement accuracy: $\pm 2\%$,
Automatic measurements: 1 - 60 min, 2 - 4 h, STAT (5 minutes continuous measurement),
SpO2:
Pulse measurement range: 20 - 300 beats per minute,
Pulse update: every 1 s,
Accuracy of saturation measurement for values 70 - 99%: $\pm 2\%$,
Equipped with the following modules and cabling:
SPO2 - Clip sensor,
NIBP - Unimed cuff size 25 - 35 cm with Spacelabs air hose,
9-pin male COM serial port connector,
Patient monitor placed on a mobile trolley:
Steel monitor stand with secure screw fixation,
Rubberised carrying handle,
Five-star base on Ø 8 cm castors, two of which with brakes,
Included:
English user manual (PDF),
Power cord,
Dimensions:
Monitor: 19 x 23 x 17 cm,
Monitor on stand: Ø 62 x 134 cm,
Weight: 7.5 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

