



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



INJECTA sterilisation box (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 1500 W,

INJECTA sterilization box is designed for dry sterilization of instruments and their subsequent storage in sterile conditions,

The can consists of:

A tub with a heater,

Perforated "sieve" insert with handles,

Sealed lid with handles,

Product features:

Has accidental opening protection in the form of two overlapping handles,

Made of steel,

Ceramic handles and feet,

Easy to keep clean,

Includes:

2 handles for carrying the insert,

Power cord,

Dimensions:

Box: 19 x 18 x 54 cm,

Perforated insert: 18 x 4 x 42 cm,

Weight: 5.3 kg,

It has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

