



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Olympus Polypectomy Loop SD - 6L AE (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Japanese manufacture,  
Olympus SD - 6 is a hexagonal loop with handle for electrosurgical removal of polyps under endoscopic observation,  
Product features:  
3 - ring ergonomic handle grip guarantees the operator a firm grip, control and maneuverability,  
Easily grasps polyps even in hard to reach areas,  
Loops ensure precise cutting,  
Loop requires power supply: Model PSD diathermy or Model UES electrosurgical unit (PSD-3, UES-2, UES-10, PSD-2, PSD-2E, PSD-10),  
The loop consists of:  
4 hexagonal loops with handle,  
Original wire A,  
4 wire guides,  
4 leaders (wire leader),  
2 attachments  
for cleaning,  
User manual in English (PDF),  
Brochure "Radio frequency cutting and coagulation in endoscopic use; Explanation",  
Overall dimensions of the handle: 6.5 x 4 x 17 cm,  
Weight of the set: 0.3 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

