



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Gynaecological chair TEYCO MED FRANCY NEW (New)**

New,  
Made in Italy,  
Features and benefits of the armchair:  
Elegantly and comfortably designed armchair with a wide range of colours for individual configuration.  
Professional chair used in gynaecology and urology,  
Quiet and smooth adjustment of chair functions.  
Chair equipped with 4 low voltage electric motors (24V-DC),  
Frame made of high-grade steel, powder-coated in light grey with antibacterial properties,  
Fitted with four castors with brakes,  
Elements of the bed are covered with upholstery with antibacterial, antifungal properties; the material is fireproof, disinfectant and highly abrasion resistant,  
Minimum height allows comfortable access to the chair without the use of a stepladder; this includes elderly patients or patients with hip prosthesis,  
The management of the electric movements is organised by both a hand control and a foot control. Both use wireless technology.  
Standard equipment:  
Stainless steel bowl for collecting liquids: equipped with a support with articulated arm,  
Several side holders for accessories on the seat,  
Remote control with position memory function (allows 4 positions to be memorised):  
2 positions are already permanently set: CPR position and gynaecological examination starting position,  
2 positions can be freely set by the doctor (foot) managed by a switch with wireless technology. Management of electric movements is organised by a foot switch with wireless technology,  
As additional options are available (for an extra charge):  
A pair of padded knee pads made of polyurethane foam. Adjustable and equipped with articulated support with locking mechanism (12500005/A),  
Detachable ultrasound footrest (12900011),  
Paper roll holder (12900017),  
Adjustable headrest (12500007),  
Handy instrument tray (12900018),  
Rechargeable battery (12500009),  
Drip holder with 2 hooks (12000012),  
Colposcopic support integrated into the seat (60000001),  
Technical data:  
Supply voltage: 230V / 50Hz,  
Motor supply voltage: 24 V DC - IP54,  
Adjustment of seat elements:  
Height adjustment: electric,  
Backrest adjustment: electric  
Trendelenburg / Anti Trendelenburg: electrically,  
Electric control management: Electric controls are managed via a control panel and pedals with wireless technology,  
Maximum load capacity: 230 kg,  
Electrical protection factor: IP54,  
Seat weight: 135 kg,  
Minimum height of the seat (from the ground): 510 mm,  
Maximum seat height (from the ground): 960 mm,  
Holds a valid Technical Passport,  
Warranty: 24 months,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

