



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Sterilization cassette for the STATIM autoclave (None)

New, after sales,
Canadian production,
Sterilization cassette for STATIM autoclave fits 5000, 5000S models,
The cassette contains a gasket and a wire frame,

Weight: 7.2 kg,

Dimensions of cassette:

Width: 180 mm,

Length: 370 mm,

Depth: 65 mm,

It has a valid Technical Passport,

Warranty: 6 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

