



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ADIVIVE closed system for fat grafting (Reconditioned)**

Used,

Technical condition: very good,

Visual condition: very good,

Korean production, (2014),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 450 VA,

Adivive enables the collection and administration of fat tissue in a closed, sterile procedure, which has found application in procedures such as body and facial modelling, facial rejuvenation, scar filling,

The material for autotransplantation is taken with a special syringe from the areas of knees, abdomen, thighs, buttocks, and then the most valuable part is separated using the system,

The final product of Adivive is highly concentrated transplant material: young fat cells and stem cells or preadipocytes,

Main advantages of the system:

Maximum safety of treatments,

High survival rate of the transplanted tissue,

Fast and simple procedure,

Natural effect of the treatment,

Improved skin quality,

Alternative to implants,

System features:

Provides the generation of high centrifugal force during an 8 - 10 minute centrifugation time,

Pressure on the aspirate causes stratification of the fraction into 3 layers: blood with tumescent fluid, condensed fat cells, old adipocytes damaged in the process,

Automatic RCF / RPM conversion,

Possibility to program RCF / RPM and time,

Lid lock and automatic release at the end of operation for safety,

Emergency release mechanism of the door lock in case of a power failure or interruption allows samples to be taken,

Lid drop protection protects the operator and the samples during loading and unloading,

The outer part of the door is made of ABS plastic and the inner part of durable steel,

Autoclavable and corrosion resistant rotors protect samples, operator and laboratory environment from contamination,

Imbalance detection function,

SOFT braking function for slow starting and braking,

Foiled control panel,

Keypad on control panel in English,

Compatible with 60 mL Lipokit/Adivive FPU syringes,

Allows you to obtain 25 mL of finished material from a single FPU syringe,

Specifications:

10 programs,

Programming of the centrifugation time: up to 100 minutes,

Centrifugation speed (RPM): 4000 rpm,

Relative centrifugal force (RCF):  $2848 \times g$ ,

Noise level:  $\leq 52$  dB,

Included:

Medikan Centrifuge Model 416D,

Set of accessories shown in the last pictures,

Power cord,

Dimensions: 38 x 31.5 x 48 cm,

Weight of the centrifuge: 20.5 kg,

It has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

