



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET ALPHASTAR PLUS 1132.13 Operating O/R table (Reconditioned)

- Used,
- Technical condition: perfect,
- Visual condition: perfect,
- German production,
- Power supply: 230 V and battery powered via 2 new built-in batteries: 12V / 14Ah,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- The Maquet Alphastar Plus operating table was created in close collaboration with physiotherapists, surgeons, and nursing staff. It simplifies daily work in the operating room and promotes patient safety and comfort. It can be easily modified to suit the needs of any surgical discipline in just a few quick steps.
- 5-segment,
- New Memory Visco upholstery with medical certificate, in black,
- Universal countertop,
- Electrically controlled by hydraulic drive,
- Adjustable in many planes using a wireless remote control or via a panel on the column:
 - ON/OFF,
 - LOCK/UNLOCK - locking and unlocking the table,
 - NORMAL/REV - patient orientation on the operating table,
 - TREND/REV - changing the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg position within the range of $\pm 30^\circ$,
 - UP/DOWN - operating surface height from 68.5 to 107 cm (up/down)
 - TILT L/TILT R - lateral tilt in the range of 18° ,
 - BACK UP/BACK DN - adjustment of the back section position,
 - LEG UP/LEG DN - adjusting the leg segment: up/down,
 - FLEX/REFLEX - bending the tabletop in half in both directions: roof position / V position (Beach chair),
 - LEVEL - reset: return to the horizontal position of the bed,
- Manually adjustable:
 - Headrest,
 - Two-segment divided leg segment - removable, with manual side abduction,
- The base housing is made of durable plastic,
- The table has X-RAY transparency,
- Equipped with rails for mounting additional equipment,
- Table adapted to C-arm,
- Mobile, on double wheels,
- Included in the set:
 - Wireless remote control,
 - Remote control charger,
 - Power cord,
- The table can be equipped with other accessories from our offer for an additional fee, examples are shown in the last photos,
- Bed dimensions: 52 x 232 x 6 cm,
- Overall dimensions: 59 x 234 x 76/114 cm,
- Maximum load capacity: 450 kg,
- Weight: 336 kg,
- Has a valid Technical Passport issued,
- Guarantee:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Financing options (Poland only): Installments, Leasing, Loan,

DESCRIPTION AND ILLUSTRATIVE PHOTOS OF THE TABLE

