



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Set of FIAB clip-on paediatric electrodes for ECG (New)

New,

Technical condition: very good,

Visual condition: very good,

Made in Italy, (2016),

Fiab electrodes are used for ECG examinations in children,

The sensors of the electrodes made of silver chloride, which allows a clean signal transmission,

The electrodes have a connector for banana and napa type plugs,

Included:

4 pcs pediatric clip-on limb electrodes,

User manual in languages: EN, DE, IT, FR, ES, PT, RU,

Dimensions: 9 x 2 x 5 cm,

Weight: 0.1 kg,

It has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

